



# AEC-NASA TECH BRIEF



AEC-NASA Tech Briefs describe innovations resulting from the research and development program of the U.S. AEC or from AEC-NASA interagency efforts. They are issued to encourage commercial application. Tech Briefs are published by NASA and may be purchased, at 15 cents each, from the Clearinghouse for Federal Scientific and Technical Information, Springfield, Virginia 22151.

## Ceric and Ferrous Dosimeters Show Precision for 50–5000 Rad Range

### The problem:

To develop ferrous and ceric dosimeters for use in mammalian radiobiology that would be sensitive in the biologically interesting 50 to 1500 rad range. The classical ferrous sulfate (Fricke) dosimeters have been used to date; however, they exhibit low sensitivity in the desired range, have a tendency to autooxidation, and require a temperature-controlled ultraviolet spectrophotometer. In developing the improved dosimeter, systems involving the addition of sensitizers (e.g., benzoic acid) before irradiation were avoided since this created a different radiochemical situation and required redetermination of radiochemical yields (G values) as functions of linear energy transfer (LET), dose rate, etc. Similarly, determinations of remaining ferrous ions were avoided, for they led to poor precision and reproducibility, especially if exposure lasted a period of several hours.

### The solution:

Ammonium thiocyanate, added to the usual ferrous sulfate dosimeter solution, yielded a very stable, precise and temperature-independent system—eight times as sensitive as the classical Fricke system in the 50 to 5000 rad range.

Ceric dosimeters, which are promising for use in mixed radiation fields, possessed a response nearly independent of LET. Fluorometry of an irradiated ceric solution, to which fluorescein was then added, afforded a method several hundred times more sensitive than previous ceric dosimeters.

### How it's done:

The ferrous dosimeter was produced by adding ammonium thiocyanate to a ferrous sulfate dosimeter

solution which had been irradiated and handled in the usual way. It was composed of 0.001 M in  $\text{Fe}(\text{NH}_4)_2(\text{SO}_4)_2$ , 0.001 M in NaCl, 0.4 M in  $\text{H}_2\text{SO}_4$ , and was aerated before use.

This combination produced a two-fold increase in sensitivity over previous thiocyanate methods because of the increased SCN concentration obtainable with the ammonium salt  $\text{NH}_4\text{SCN}$ , also diminished temperature sensitivity, and decreased the rate of ferrous ion autooxidation.

Two methods have been used successfully for the determination of the ferric thiocyanate complex. Both employ visible region colorimeters at 4750 and 4820 Å. In the first, one volume of a saturated aqueous solution of  $\text{NH}_4\text{SCN}$  was added to three volumes of irradiated dosimeter. In the second, dosimeter was added to a known mass of solid ammonium thiocyanate. This latter method was more sensitive because of the lack of dilution.

Absorbance was strictly proportional to ferric ion concentration, with  $\epsilon_{4820} = 17,600 \text{ liter mole}^{-1} \text{cm}^{-1}$ , from  $A=0$  to  $A=1.0$ . No temperature effect was found between 23.5° and 31° C and no time effect between 5 to 25 minutes after mixing. In a 10-cm cuvette, 50 rads gave  $A=0.126 \pm 0.001$ .

The ceric dosimeter system was based on the observation that the fluorescence of fluorescein is quantitatively destroyed by the ceric ion. Fluorescence of one volume of aqueous fluorescein, added to 12.5 volumes of irradiated ceric dosimeter, was found to be a quasi-linear function of dose over the 50–5000 rad range with a precision of  $\pm 1\%$  and no observable temperature or time dependence between 20° to 30° C and 5 to 100 minutes.

(continued overleaf)

**Notes:**

1. Reference: For additional details see: "Solid State and Chemical Radiation Dosimetry in Medicine and Biology," IAEA, Vienna, 1967, pp. 257-266. "Ceric and Ferrous Dosimeters for the 50-5000-rad Range" by N. A. Frigerio and V. D. Henry, Division of Biological and Medical Research, Argonne National Laboratory.
2. Inquiries concerning this innovation may be directed to:

Office of Industrial Cooperation  
Argonne National Laboratory  
9700 South Cass Avenue  
Argonne, Illinois 60439  
Reference: B68-10426

Source: N. A. Frigerio and V. D. Henry  
of the Biological and Medical  
Research Division  
(ARG-10173)

**Patent status:**

Inquiries about obtaining rights for commercial use of this innovation may be made to:

Mr. George H. Lee, Chief  
Chicago Patent Group  
U.S. Atomic Energy Commission  
Chicago Operations Office  
9800 South Cass Avenue  
Argonne, Illinois 60439